

HITACHI 717

APOLIPOPROTEIN A1 FS

Chemistry settings

Temperature: 37 °C

PROGRAM 2: CHEMISTRY PARAMETERS

TEST	APOA
ASSAY CODE	2 - 24 - 50
SAMPLE VOLUME (µl)	2
R 1 VOLUME (µl)	250 - 100 - NO
R 2 VOLUME (µl)	50 - 20 - NO
WAVELENGTH (nm)	0 - 570
CALIBR. METHOD	3 - 1 - 4
STD.1 CONC.-POS.	0 - 1
STD.2 CONC.-POS.	* - *
STD.3 CONC.-POS.	* - *
STD.4 CONC.-POS.	* - *
STD.5 CONC.-POS.	0 - 0
STD.6 CONC.-POS.	0 - 0
SD LIMIT	999
DUPLICATE LIMIT	500
SENSITIVITY LIMIT	0
ABS.LIMIT(INC/DEC)	0 - 0
PROZONE LIMIT	32000 - 1
EXPECTED VALUE	110 - 190
PANIC VALUE	0.2 - 250
INSTRUMENT FACTOR	1.0

#) Data entry by the user

*) Enter calibration or standard value and position

Order information

Cat. No.	Kit size
1 7102 99 10 015	R1 2 x 25 ml + R2 1 x 10 ml
10 710 021	R1 5 x 25 ml + R2 1 x 25 ml
1 7102 99 10 730	R1 4 x 20 ml + R2 2 x 8 ml
1 7100 99 10 041	3 x 1 ml TruCal Apo A1: Calibrator set with 3 different levels

Notes

1. Please refer to the package insert for Apolipoprotein A1 FS for detailed information about the test on the following:

- Clinical Relevance
- Method and Principle
- Composition and Stability of the Reagents
- Specimens
- Calibrators and Controls
- Performance Characteristics regarding
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
 - Method Comparison
- Reference Ranges
- Literature

2. The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
3. Manufactured by
DiaSys Diagnostic Systems GmbH & Co.KG
Alte Strasse 9, 65558 Holzheim, Germany