

BM / HITACHI 912

04-01 CHEMISTRY PARAMETERS

TEST [APOA] [#]	TEST NAME [APOA]	UNIT [mg/dl]							
DATA MODE [ON BOARD]	REPORT NAME [APO A1]								
CONTROL INT [0]	INSTR. FACT. (Y=aX + b) a [1.0]	b [0.0]							
EXPECTED VALUE	CLASS 1	EXPECTED VALUE CLASS 2							
AGE	M	F							
[] [] []	- [] []	- [] []							
[] [] []	- [] []	- [] []							
	[110] - [170]	[120] - [190] []							
TECHNICAL LIMIT	CLASS I	CLASS 2							
	[0.2] - [250]	[] - []							
STD	CON	POS	S.VOL	PRE.	DIL.	VOL	CODE	LOT	QUALITATIVE
									[NO]
(1)	[0.0]	[#]	[2]	[0]	[0]	[#]			(1) [] []
(2)	[*]	[#]	[2]	[0]	[0]	[#]			(2) [] []
(3)	[*]	[#]	[2]	[0]	[0]	[#]			(3) [] []
(4)	[*]	[#]	[2]	[0]	[0]	[#]			(4) [] []
(5)	[]	[#]	[]	[0]	[0]	[#]			(5) [] []
(6)	[]	[#]	[]	[0]	[0]	[#]			(6) [] []

TEST [APOA]
ASSAY CODE [2 point end] [10] [] WAVELENGTH (SUB / MAIN) [700] / [580]
ASSAY POINTS [15]- [31] - [0]- [0] DILUTION [301] [99] < CLASS 1 > < CLASS 2 >
S.VOL (NORMAL) [2] [0] [0] [] [] []
S.VOL (DECREASE) [1] [0] [0] [] [] []
S.VOL (INCREASE) [4] [0] [0] [] [] []
ABS. LIMIT [0] [] [2:INCREASE]
PROZONE LIMIT [32000] [] [1:HIGHER]
REAGENT R1 [250] [0] [#] [#]
R2 [0] [0] [#] [#]
R3 [50] [0] [#] [#]
R4 [0] [0] [#] [#]
CALIB. TYPE [NON-LINEAR] [4] [4] [0] []
AUTOCALIB.
TIME OUT BLANK [0] SD LIMIT [999]
SPAN [0] DUPLICATE LIMIT [500]
2 POINT [0] SENSITIVITY LIMIT [0]
FULL [0] SI ABS. LIMIT [-32000][32000]
CHANGE LOT [NO] COMPENSATED LIMIT []
BOTTLE [NO]

Enter the next code

APOLIPOPROTEIN A1 FS

Order information

Cat. No.	Kit size
1 7102 99 10 015	R1 2 x 25 ml + R2 1 x 10 ml
10 710 021	R1 5 x 25 ml + R2 1 x 25 ml
1 7102 99 10 730	R1 4x 20 ml + R2 2x 8 ml
1 7100 99 10 041	3 x 1 ml TruCal Apo A1: Calibrator set with 3 different levels

Notes

- Please refer to the package insert for Apolipoprotein A1 FS for the detailed information about the test on the following:

- Clinical Relevance
- Method and Principle
- Composition and Stability of the Reagents
- Specimens
- Calibrators and Controls
- Performance Characteristics concerning:
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
 - Method Comparison
- Reference Ranges
- Literature

- The stability of the reagent on board of the analyser is at least one month provided that contamination and evaporation are avoided.

- Manufactured by
DiaSys Diagnostic Systems GmbH & Co.KG
Alte Strasse 9, 65558 Holzheim, Germany.

Data entry by the user
 * Enter calibration or standard value
 ** The given factor must be checked by a calibration serum.