

# COBAS MIRA Plus

# APOLIPOPROTEIN B FS

## Instrument settings

Test name: ApoB

Temperature: 37 °C

<b>GENERAL</b>			
Measurement mode:	ABSORB		[1]
Reaction mode:	R - S - SR1		[3]
Calibration mode:	LOGIT / LOG 4		[5]
Reagent blank:	REAG / DIL		[2]
Cleaner:	NO		[1]
Wavelength:	340 nm		[1]
Decimal position:	0		
Unit:			[12 (mg / dl)]
<b>ANALYSIS</b>			
Post dil. factor:			[2]
conc. factor:	NO		[SPACE]
Sample cycle:			[1]
volume:			3 µl
Dilution name:	H2O		[00]
volume:			10 µl
Reagent cycle:			1
volume:			250 µl
Start R 1 cycle:			6
volume:			50 µl
Dilution name:			(00)
volume:			20 µl
<b>CALCULATION</b>			
Sample limit:	NO		[SPACE]
Point:			
Reac. direction:	INCREASE		[1]
Check:	ON		[1]
Convers. factor:			1
Offset:			0
Test range low:		or	25 (mg/dl)
high:		or	250 (mg/dl)
Norm. range low:		or	45 (mg/dl)
high:		or	135 (mg/dl)
Number of steps:			1
Calc. step A:	ENDPOINT		[1]
Readings first:			5
last:			19
Reaction limit:			
Point:			
<b>CALIBRATION</b>			
Calibration interval:	ON REQUEST		[3]
Time:	NO		[SPACE]
Blank:			
Reag. range low:	NO		[SPACE]
high:	NO		[SPACE]
Blank range low:	NO		[SPACE]
high:	NO		[SPACE]
Factor:			
Standard pos:			# ....
std.-1:	* Calibrator value		* ....
std.-2:	* Calibrator value		[* ....]
std.-3:	* Calibrator value		[* ....]
Replicate:	TRIPL		[3]
Deviation:			[5 %]
<b>CONTROL</b>			
CS 1 Pos:	Low:.....	Assign:....	High:.....
CS 2 Pos:	Low:.....	Assign:....	High:.....
CS 3 Pos:	Low:.....	Assign:....	High:.....

## Order information

Cat. No.	Kit size
1 7112 99 10 015	R1 2 x 25 ml + R2 1 x 10 ml

## Notes

- Please refer to the package insert for Apolipoprotein B FS for detailed information about the test on the following:

Clinical Relevance  
 Method and Principle  
 Composition and Stability of the Reagents  
 Specimens  
 Calibrators and Controls  
 Performance Characteristics regarding  
 - Measuring Range  
 - Specificity/Interferences  
 - Sensitivity/Limit of Detection  
 - Precision (Reproducibility, Repeatability)  
 - Method Comparison  
 Reference Ranges  
 Literature

- The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by  
 DiaSys Diagnostic Systems GmbH & Co.KG  
 Alte Strasse 9, 65558 Holzheim, Germany

- #) Data entry by the user  
 \*) Enter calibration or standard value and position  
 \*\*) The factor must be checked by a calibration serum.