

BM / HITACHI 911

04-01 CHEMISTRY PARAMETERS

TEST	[APOB]	[#]	TEST NAME	[APOB]	UNIT	[mg/dl]	
DATA MODE	[ON BOARD]		REPORT NAME	[APO B]			
CONTROL INT	[0]		INSTR. FACT.	(Y=aX + b)	a	[1.0]	
					b	[0.0]	
EXPECTED VALUE	CLASS 1		EXPECTED VALUE	CLASS 2			
AGE	M		F				
[]	[]	[]	- []	[]	- []	[]	
[]	[]	[]	- []	[]	- []	[]	
	[80]	- [155]	[75]	- [150]	[]		
TECHNICAL LIMIT	CLASS 1		CLASS 2				
	[0.3]	- [250]	[]	- []			
STD	CO	PO	S.VO	PRE. DIL.	VOL	CODE	
	N	S	L.			LOT	
						QUALITATIVE	
						[NO]	
(1)	[0.0]	[#]	[3]	[0]	[0]	[#]	(1) [] []
(2)	[*]	[#]	[3]	[0]	[0]	[#]	(2) [] []
(3)	[*]	[#]	[3]	[0]	[0]	[#]	(3) [] []
(4)	[*]	[#]	[3]	[0]	[0]	[#]	(4) [] []
(5)	[]	[#]	[]	[0]	[0]	[#]	(5) [] []
(6)	[]	[#]	[]	[0]	[0]	[#]	(6) [] []

TEST	[APOB]					
ASSAY CODE	[2 point end]	[10]	[]	WAVELGTH	(SUB / MAIN)	
						[700] / [340]
ASSAY POINTS	[15]- [31]	- [0]- [0]	DILUTION	[301]	[99]	
	< CLASS 1 >	< CLASS 2 >				
S.VOL (NORMAL)		[3]	[0]	[0]	[]	[]
S. VOL (DECREASE)		[1]	[0]	[0]	[]	[]
S.VOL (INCREASE)		[6]	[0]	[0]	[]	[]
ABS. LIMIT	[0]		[]		[2:INCREASE]	
PROZONE LIMIT	[32000]		[]		[1:HIGHER]	
REAGENT	R1	[250]	[0]	[#]	[#]	
	R2	[0]	[0]	[#]	[#]	
	R3	[50]	[0]	[#]	[#]	
	R4	[0]	[0]	[#]	[#]	
CALIB. TYPE	[NON-LINEAR]		[4]	[4]	[0]	[]
AUTOCALIB.						
TIME OUT BLANK	[0]		SD LIMIT		[999]	
SPAN	[0]		DUPLICATE LIMIT		[500]	
2 POINT	[0]		SENSITIVITY LIMIT		[0]	
FULL	[0]		SI ABS. LIMIT		[-32000][32000]	
CHANGE LOT	[NO]		COMPENSATED LIMIT		[]	
BOTTLE	[NO]					

- # Data entry by the user
- * Enter calibration or standard value
- ** The given factor must be checked by a calibration serum.
- ## Enter the next code

APOLIPOPROTEIN B FS

Order information

Cat. No.	Kit size
1 7112 99 10 015	R1 2 x 25 ml + R2 1 x 10 ml
10 711 021	R1 5 x 25 ml + R2 1 x 25 ml
1 7112 99 10 730	R1 4 x 20 ml + R2 2 x 8 ml
1 7110 99 10 041	3 x 1 ml TruCal Apo B: Calibrator set with 3 different levels

Notes

- Please refer to the package insert for Apolipoprotein B FS for detailed information about the test on the following:

- Clinical Relevance
- Method and Principle
- Composition and Stability of the Reagents
- Specimens
- Calibrators and Controls
- Performance Characteristics regarding
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
 - Method Comparison
- Reference Ranges
- Literature

- The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by
DiaSys Diagnostic Systems GmbH & Co.KG
Alte Strasse 9, 65558 Holzheim, Germany