

HITACHI 717

BILIRUBIN AUTO TOTAL FS

Chemistry settings

Temperature: 37 °C

PROGRAM 2: CHEMISTRY PARAMETERS

TEST	BILI
ASSAY CODE	2 - 24 - 50
SAMPLE VOLUME (µl)	6
R 1 VOLUME (µl)	240 - 100 - 0
R 2 VOLUME (µl)	60 - 20 - 0
WAVELENGTH (nm)	660 - 546
CALIBR. METHOD	1 - 0 - 0
STD.1 CONC.-POS.	0.0 - 1
STD.2 CONC.-POS.	* - *
STD.3 CONC.-POS.	0 - 0
STD.4 CONC.-POS.	0 - 0
STD.5 CONC.-POS.	0 - 0
STD.6 CONC.-POS.	0 - 0
SD LIMIT	0.1
DUPLICATE LIMIT	100
SENSITIVITY LIMIT	0
ABS.LIMIT(INC/DEC)	0 - 0
PROZONE LIMIT	0 - 0
EXPECTED VALUE	0.1 - 1.2
PANIC VALUE	0.1 - 30.0
INSTRUMENT FACTOR	1.0

#) Data entry by the user

*) Enter calibration or standard value and position

Order information

Cat. No.	Kit size					
10 081 021	R1	5 x	20 ml	+	R2	1 x 25 ml
10 081 022	R1	5 x	80 ml	+	R2	1 x 100 ml
10 081 023	R1	1 x	800 ml	+	R2	1 x 200 ml
10 081 717	R1	5 x	80 ml	+	R2	5 x 20 ml

Notes

1. Please refer to the package insert for Bilirubin Auto Total FS for detailed information about the test on the following:

Clinical Relevance
Method and Principle
Composition and Stability of the Reagents
Specimens
Calibrators and Controls
Performance Characteristics regarding
- Measuring Range
- Specificity/Interferences
- Sensitivity/Limit of Detection
- Precision (Reproducibility, Repeatability)
- Method Comparison
Reference Ranges
Literature

2. The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
3. Manufactured by
DiaSys Diagnostic Systems GmbH & Co.KG
Alte Strasse 9, 65558 Holzheim, Germany