

BILIRUBIN AUTO TOTAL FS*

Order information

Cat. No.	Kit size								
10 081 021	R1 5 x	20 ml	+	R2 1 x	25 ml				
10 081 022	R1 5 x	80 ml	+	R2 1 x	100 ml				
10 081 023	R1 1 x	800 ml	+	R2 1 x	200 ml				

Method

Colorimetric test, "DCA (2,4-Dichloroaniline)"

Total Bilirubin in presence of diazotized 2,4-Dichloroaniline forms a red colored azocompound in acidic solution. A specific mixture of detergents enables a safe determination of the Total Bilirubin.

Reagent preparation and stability

The reagents are ready-to-use and stable up to the end of the indicated month of expiry, if contamination is avoided and stored at 2 – 8 °C. The reagent 2 must be protected from light.

Specimen

Serum, heparinized or EDTA plasma. Avoid hemolysis!
Store protected from light.

Stability:	1 day	at	15 – 25 °C
	4 days	at	2 – 8 °C
	up to 3 months	at	- 20 °C
	(not in case of repeated deep freezing)		

Components and concentration in the test

R1:	TRIS	pH 8.2	8 mmol/l
	NaCl		7 g/l
	Detergents		
R2:	2,4-Dichlorophenyl Diazoniumsalt		1 mmol/l
	HCl		30 mmol/l
	Detergents		

Notes

Ascorbic acid and Hemoglobin interfere.

Normal range (see reference 2,3)

Newborn	0 - 24 h:	up to	5 mg/dl (86 µmol/l)
	24 - 48 h:	up to	9 mg/dl (155 µmol/l)
	3 – 5 days:	up to	12 mg/dl (205 µmol/l)
Adults:	after 4 weeks:	up to	1.5 mg/dl (26 µmol/l)
		up to	1.1 mg/dl (18.8 µmol/l)

References

- Rand, R. N., di Pasqua, A., Clin Chem., 8, (1962), 570
- Weigl, E., Bach, H., Krieg, D., Med. Klin., 70, (1975), 664 – 669
- Keller, H.: Klinisch-chemische Labordiagnostik für die Praxis, 2nd edition, Georg Thieme Verlag, Stuttgart 1991, 246

HITACHI 917

Chemistry parameters

Analysis		TBIL				Ser/Pl			
Test / Type		2 Point End	A	10	A	15	34	0	0
Assay / Time / Point		660	A	546	A				
Wave (2 nd /Primary)		5	0	0					
S.Vol (Normal)		1	0	0					
S.Vol (Decrease)		10	0	0					
S.Vol (Increase)		00951	99						
Diluent		200	0	#	60				Timing
Reagent (R1) T1		0	0	#	0				R1
Reagent (R2) T2		50	0	#	60				R3
Reagent (R3) T3		0	0	#	0				
Reagent (R4) T4		12000	Increase		A				
Abs. Limit		32000	0	Lower	A				
Prozone Limit		Detergent 1		A					
Cell Detergent									

Calibration	
Calibration type	Linear A A
Point	2 Span Point 2
Weight	0
Autocalibration	
Blank	Blank A
Span	Blank A
2Point	
Full	
SD Limit	0.1
Duplicate limit	10 % 200 Abs
Sensitivity limit	-99999 999999
S1 Abs limit	-32000 32000

Range	
Application Code	# Unit mg/dl A
Report Name	Bilirubin Total
Data Mode	On Board A
Control Interval	1000
Instrument Factor (Y=aX+b)	a= 1.0 b= 0.0
Technical Limit	0.0 30.0
Repeat Limit	0.0 30.0
Expected Value	
(Male)	Y A 0.0 1.1
(Female)	Y A
(Default)	Male A Range3 A
Qualitative	Cancel A
	(1) 0
	(2) 0
	(3) 0
	(4) 0
	(5) 0
	(6) 0

Others	
<Standard>	(1) (2) (3) (4) (5) (6)
Calib. Code	501 #
Concentration	0.0 *
Position	
Sample Volume	5 5
Diluent S.Vol	0 0
Diluent Volume	0 0

#) Data entry by the user

*) Enter calibration or standard value

* fluid stable