

OLYMPUS AU 800
Individual test parameters

Bilirubin Auto Total FS

1.	TEST #	SERUM	
2.	SAMPLE VOLUME		(6)
3.	SAMPLE VOLUME for Repeat Run		(#)
4.	REAGENT VOLUME	R1	(250)
	Diluent Volume		()
	Diluent Volume	R2	(60)
			()
5.	Wavelength	1	(540)
		2	(660)
6.	Method		(END)
7.	Reaction Slope		(+)
8.	Measuring Point	S1	(0)
		E1	(16)
		S2	(0)
		E2	(6)
9.	QD Value Range		
		Max	(+2.00)
		Min	(-2.00)
10.	Limit of Linearity		(20)
11.	No-Lag-Time		(NO)
12.	Data Check 1		
		Px	(0)
		Py	(0)
		Pz	(0)
13.	Data Check 2		
		Px	(0)
		Py	(0)
		Pz	(0)
14.	Normal Value	H	(17)
		L	(0)
15.	Dynamic Range	H	(513)
		L	(0)
16.	Panic Value	H	(#)
		L	(#)
17.	Reagent QD Range		
	First Point	H	(2.00)
		L	(-2.00)
	Last Point	H	(2.00)
		L	(-2.00)
18.	Decision Levels for Repeat Run	H	(#)
		L	(#)

Order information

Cat. No.	Kit size					
10 081 021	R1 5 x	20 ml	+	R2 1 x	25 ml	
10 081 022	R1 5 x	80 ml	+	R2 1 x	100 ml	
10 081 023	R1 1 x	800 ml	+	R2 1 x	200 ml	

Notes

1. Please refer to the package insert for Bilirubin Auto Total FS for detailed information about the test on the following:

Clinical Relevance
 Method and Principle
 Composition and Stability of the Reagents
 Specimens
 Calibrators and Controls
 Performance Characteristics regarding
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
 - Method Comparison
 Reference Ranges
 Literature

2. The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
3. Manufactured by
 DiaSys Diagnostic Systems GmbH & Co.KG
 Alte Strasse 9, 65558 Holzheim, Germany

#) Data entry by the user
 *) Calculated by the analyzer