

FALCOR

BILIRUBIN AUTO TOTAL FS

DESCRIPTION / MODIFICATION OF METHOD

Type	: End Pt.
Unit	: mg/dl
Number of Decimals	: 2
Diluent	: 99 Diluan
Pre-washing	: No
First reagent Vol. (µl)	: 200
Second reagent Vol. (µl)	: 50
Sample Volume (µl)	: 10.0
Lower normal value	: 0.0
Higher normal value	: 1.1
Dilution factor	: -
Systematic dilution	: 1
Dilution for Urine	: 1
Low reagent blank	: 0.0300
High reagent blank	: 0.1000
Maximum Deviation	: 0.0200
Primary Wavelength	: 550
High Wavelength	: -
Low Wavelength	: -
	:
Depletion limit	: -
T Tolerance factor	: 0.0010
C Compensation factor	: -
Slope	: 1.00
Intercept	: 0.0
Delay for reagent 2 dispense	: 02:36
Reaction delay	: 03:24
Incubation	: -
Total	: -

Order information

Cat. No.	Kit size
10 081 021	R1 5 x 20 ml + R2 1 x 25 ml
10 081 022	R1 5 x 80 ml + R2 1 x 100 ml
10 081 023	R1 1 x 800 ml + R2 1 x 200 ml
10 081 704	R1 8 x 50 ml + R2 8 x 12.5 ml
10 081 717	R1 5 x 80 ml + R2 5 x 20 ml

Notes

1. Please refer to the package insert for Bilirubin Auto Total FS for detailed information about the test on the following:

Clinical Relevance
Method and Principle
Composition and Stability of the Reagents
Specimens
Calibrators and Controls
Performance Characteristics concerning
- Measuring Range
- Specificity/Interferences
- Sensitivity/Limit of Detection
- Precision (Reproducibility, Repeatability)
- Method Comparison
Reference Ranges
Literature

2. The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
3. Manufactured by
DiaSys Diagnostic Systems GmbH & Co.KG
Alte Strasse 9, 65558 Holzheim, Germany.

- #) Data entry by the user
*) Calculated by the analyzer