

MICROLAB 200

TEST PARAMETERS

1	NAME	T. BILIRUBIN
2	MODE	ENDPOINT
3	FILTER	546
4	TEMP	37°C
5	UNITS	mg/dl
6	R-BLANK	YES
	-L-Abs	0.000
	-H-Abs	2.300
1	CAL-NAME	TRUCAL
	.NR	1
	.CONC	*
	.REP	1
	.L-Abs	0.000
	.H-Abs	2.300
2	FACTOR	CALCULATED
3	S-BLANK	NO
1	DELAY	2
2	MEAS TIME	2
3	VOLUME	500
1	TEST REP	1
	.MAX DEV	0
2	REF-LOW	0.1
3	REF-HIGH	1.2
4	L-Abs	0.000
5	H-Abs	2.300
1	CONTROL-1	#
	.BATCH-NR	#
	.TARGET	#
	.LOW-LIM	#
	.HIGH-LIM	#
2	CONTROL-2	#
	.BATCH-NR	#
	.TARGET	#
	.LOW-LIM	#
	.HIGH-LIM	#

#) DATA ENTRY BY THE USER

*) ENTER CALIBRATION OR STANDARD VALUE AND POSITION

**) THE FACTOR MUST BE CHECKED BY A CALIBRATION SERUM.

BILIRUBIN AUTO TOTAL FS

Order information

Cat. No.	Kit size
10 081 021	R1 5 x 20 ml + R2 1 x 25 ml
10 081 022	R1 5 x 80 ml + R2 1 x 100 ml
10 081 023	R1 1 x 800 ml + R2 1 x 200 ml

Notes

- Please refer to the package insert for Bilirubin Auto Total FS for the detailed information about the test on the following:

Clinical Relevance
 Method and Principle
 Composition and Stability of the Reagents
 Specimens
 Calibrators and Controls
 Performance Characteristics concerning:
 Measuring Range
 Specificity/Interferences
 Sensitivity/Limit of Detection
 Precision (Reproducibility,
 Repeatability)
 Method Comparison
 Reference Ranges
 Literature

- The stability of the reagent on board of the analyser is at least one month provided that contamination and evaporation are avoided.

- Manufactured by
 DiaSys Diagnostic Systems GmbH & Co.KG
 Alte Strasse 9, 65558 Holzheim, Germany.