

HITACHI 717

CK-MB FS

Chemistry settings

Temperature: 37 °C

PROGRAM 2: CHEMISTRY PARAMETERS

TEST	CKMB
ASSAY CODE	5 – 40 – 50
SAMPLE VOLUME (µl)	10
R 1 VOLUME (µl)	200 – 100 – NO
R 2 VOLUME (µl)	50 – 20 – NO
WAVELENGTH (nm)	405 – 340
CALIBR. METHOD	**2 – 0 – 0
STD.1 CONC.-POS.	0 – 1
STD.2 CONC.-POS.	0 – 0
STD.3 CONC.-POS.	0 – 0
STD.4 CONC.-POS.	0 – 0
STD.5 CONC.-POS.	0 – 0
STD.6 CONC.-POS.	0 – 0
SD LIMIT	0.1
DUPLICATE LIMIT	100
SENSITIVITY LIMIT	0
ABS.LIMIT(INC/DEC)	9000 – 0
PROZONE LIMIT	0
EXPECTED VALUE	0 – 24
PANIC VALUE	0 - 1000
INSTRUMENT FACTOR	1.0

#) Data entry by the user

*) Enter calibration or standard value and position

**) Factor must be checked by a calibration serum

Order information

Cat. No.	Kit size
10 165 021	R1 5 x 20 ml + R2 1 x 25 ml
10 165 022	R1 5 x 80 ml + R2 1 x 100 ml
10 165 023	R1 1 x 800 ml + R2 1 x 200 ml
10 165 700	R1 16 x 20 ml + R2 4 x 20 ml

Notes

1. Please refer to the package insert for CK-MB FS for detailed information about the test on the following:

Clinical Relevance
Method and Principle
Composition and Stability of the Reagents
Specimens
Calibrators and Controls
Performance Characteristics regarding
- Measuring Range
- Specificity/Interferences
- Sensitivity/Limit of Detection
- Precision (Reproducibility, Repeatability)
- Method Comparison
Reference Ranges
Literature

2. The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
3. Manufactured by
DiaSys Diagnostic Systems GmbH & Co.KG
Alte Strasse 9, 65558 Holzheim, Germany