

HITACHI 902

INSTRUMENT SETTINGS

CK-MB FS

No.	<Chemistry>	
1	Test Name	CKMB
2	Assay Code (Mthd)	Rate A
3	Assay Code (2. Test)	0
4	Reaction Time	10
5	Assay Point 1	28
6	Assay Point 2	35
7	Assay Point 3	0
8	Assay Point 4	0
9	Wave Leng. (SUB)	415
10	Wave Leng. (MAIN)	340
11	Sample Volume	10
12	R1 VOLUME	200
13	R1 Pos.	#
14	R1 Bottle Size	Large
15	R2 VOLUME	0
16	R2 Pos.	0
17	R2 Bottle Size	Small
18	R3 VOLUME	50
19	R3 Pos.	#
20	R3 Bottle Size	Small
21	Calib. Type (Type)	Linear
22	Calib. Type (Wght)	0
23	Calib. Conc. 1	0
24	Calib. Pos. 1	#
25	Calib. Conc. 2	#
26	Calib. Pos. 2	#
27	Calib. Conc. 3	0
28	Calib. Pos. 3	0
29	Calib. Conc. 4	0
30	Calib. Pos. 4	0
31	Calib. Conc. 5	0
32	Calib. Pos. 5	0
33	Calib. Conc. 6	0
34	Calib. Pos. 6	0
35	S 1 ABS.	0
36	K Factor	10000
37	K 2 Factor	10000
38	K 3 Factor	10000
39	K 4 Factor	10000
40	K 5 Factor	10000
41	A Factor	0
42	B Factor	0
43	C Factor	0
44	SD Limit	0.1
45	Duplicate Limit	300
46	Sens. Limit	0
47	S 1 ABS Limit (L)	-32000
48	S 1 ABS Limit (H)	32000
49	ABS Limit	9000
50	ABS Limit (D/I)	Increase
51	Prz. Limit	0
52	Prz. Limit (U/D)	Lower
53	Prz. (End Point)	35
54	Expect. Value (L)	0
55	Expect. Value (H)	24
56	Instr. Fact. (a)	1
57	Instr. Fact. (b)	0
58	Key Setting	#

Order information

Cat. No. Kit size
10 165 021 R1 5 x 20 ml + R2 1 x 25 ml

Notes

1. Please refer to the package insert for CK-MB FS for the detailed information about the test on the following:

Clinical Relevance
Method and Principle
Composition and Stability of the Reagents
Specimens
Calibrators and Controls
Performance Characteristics concerning:
Measuring Range
Specificity/Interferences
Sensitivity/Limit of Detection
Precision (Reproducibility,
Repeatability)
Method Comparison
Reference Ranges
Literature

2. The stability of the reagent on board of the analyser is at least one month provided that contamination and evaporation are avoided.

3. Manufactured by
DiaSys Diagnostic Systems GmbH & Co.KG
Alte Strasse 9, 65558 Holzheim, Germany.

#) Data entry by the user

*) Calculated by the analyzer