

BM / HITACHI 911

04-01 CHEMISTRY PARAMETERS

CK-MB FS

TEST [CKMB] [#]	TEST NAME [CKMB]	UNIT [U/l]					
DATA MODE [ON BOARD]	REPORT NAME [CK-MB]						
CONTROL INT [0]	INSTR. FACT. (Y=aX + b) a [1.0]	b [0.0]					
EXPECTED VALUE	CLASS 1	EXPECTED VALUE CLASS 2					
AGE	M	F					
[] [] []	- [] []	- [] []					
[] [] []	- [] []	[] []					
[0]	- [24]	[0] - [24] []					
TECHNICAL LIMIT	CLASS 1	CLASS 2					
[0]	- [1000]	[] - []					
STD	CO	PO	S.VO	PRE. DIL.	VOL	CODE	QUALITATIVE
	N	S	L.			LOT	[NO]
(1)	[0.0]	[#]	[10]	[0]	[0]	[#]	(1) [] []
(2)	[*]	[#]	[10]	[0]	[0]	[#]	(2) [] []
(3)	[]	[]	[]	[]	[]	[]	(3) [] []
(4)	[]	[]	[]	[]	[]	[]	(4) [] []
(5)	[]	[]	[]	[]	[]	[]	(5) [] []
(6)	[]	[]	[]	[]	[]	[]	(6) [] []

Order information

Cat. No.	Kit size						
10 165 021	R1 5 x	20 ml +	R2 1 x	25 ml			
10 165 022	R1 5 x	80 ml +	R2 1 x	100 ml			
10 165 023	R1 1 x	800 ml +	R2 1 x	200 ml			
10 165 700	R1 16 x	20 ml +	R2 4 x	20 ml			

Notes

- Please refer to the package insert for CK-MB FS for detailed information about the test on the following:

Clinical Relevance
 Method and Principle
 Composition and Stability of the Reagents
 Specimens
 Calibrators and Controls
 Performance Characteristics regarding
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
 - Method Comparison
 Reference Ranges
 Literature

- The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by
 DiaSys Diagnostic Systems GmbH & Co.KG
 Alte Strasse 9, 65558 Holzheim, Germany

TEST [CKMB]	
ASSAY CODE [Rate A] [10] []	WAVELGTH (SUB / MAIN) [415] / [340]
ASSAY POINTS [25] - [31] - [0] - [0]	DILUTION [301] [99]
< CLASS 1 >	< CLASS 2 >
S.VOL (NORMAL)	[10] [0] [0] [] [] []
S. VOL (DECREASE)	[5] [0] [0] [] [] []
S.VOL (INCREASE)	[20] [0] [0] [] [] []
ABS. LIMIT [9000]	[] [2:INCREASE]
PROZONE LIMIT [0]	[] [LOWER]
REAGENT R1 [200] [0] [#] [#]	
R2 [0] [0] [#] [#]	
R3 [50] [0] [#] [#]	
R4 [0] [0] [#] [#]	
CALIB. TYPE [1:LINEAR]	[2] [2] [0] []
AUTOCALIB.	
TIME OUT BLANK [0]	SD LIMIT [0.1]
SPAN [0]	DUPLICATE LIMIT [200]
2 POINT [0]	SENSITIVITY LIMIT [0]
FULL [0]	SI ABS. LIMIT [-32000][32000]
CHANGE LOT [NO]	COMPENSATED LIMIT []
BOTTLE [NO]	

- # Data entry by the user
 * Enter calibration or standard value
 ** The given factor must be checked by a calibration serum.
 ## Enter the next code