

COBAS MIRA Plus

Monoreagent Application

CK-MB FS

Instrument settings

Test name: CK-MB

Temperature: 37 °C

GENERAL		
MEASUREMENT MODE:	ABSORB	[1]
REACTION MODE:	R -S	[1]
CALIBRATION MODE:	FACTOR	[1]
REAGENT BLANK:	REAG / DIL	[1]
CLEANER:	NO	[1]
WAVELENGTH:	340 nm	[1]
DECIMAL POSITION:	0	
UNIT:	U / L	[]
ANALYSIS		
POST DIL. FACTOR:		[3]
CONC. FACTOR:	NO	[SPACE]
SAMPLE CYCLE:		[1]
VOLUME:		8 µl
DILUTION NAME:	H ₂ O	[00]
VOLUME:		20 µl
REAGENT CYCLE:		1
VOLUME:		200 µl
START R 1 CYCLE:		
VOLUME:		
DILUTION NAME:		
VOLUME:		
CALCULATION		
SAMPLE LIMIT:	0.4500	
POINT:	T1	
REAC. DIRECTION:	INCREASE	[1]
CHECK:	ON	[1]
CONVERS. FACTOR:		
OFFSET:		
TEST RANGE LOW:	0 (U / L)	
HIGH:	100 (U / L) or	
NORM. RANGE LOW:	0 (U / L)	
HIGH:	24 (U / L)	
NUMBER OF STEPS:	1	
CALC. STEP A:	Kinetic	[]
READINGS FIRST:	5	
LAST:	16	
REACTION LIMIT:	0.65	
POINT:	T1	
CALIBRATION		
CALIBRATION INTERVAL	ON REQUEST	[3]
TIME:	NO	[SPACE]
BLANK:		
REAG. RANGE LOW:	0.0500	
HIGH:	0.400	
BLANK RANGE LOW:	-0.003	
HIGH:	0.003	
FACTOR:	**14800	
STANDARD POS:		#
STD.-1:	NO	[SPACE]
STD.-2:	NO	[SPACE]
STD.-3:	NO	[SPACE]
REPLICATE:	NO	[]
DEVIATION:	NO	[]
CONTROL		
CS 1 POS:	Low:..... Assign:..... High:.....	
CS 2 POS:	Low:..... Assign:..... High:.....	
CS 3 POS:	Low:..... Assign:..... High:.....	

Order information

Cat. No.	Kit size
10 165 021	R1 5 x 20 ml + R2 1 x 25 ml
10 165 022	R1 5 x 80 ml + R2 1 x 100 ml
10 165 023	R1 1 x 800 ml + R2 1 x 200 ml

Notes

- Please refer to the package insert for CK-MB FS for detailed information about the test on the following:

Clinical Relevance
 Method and Principle
 Composition and Stability of the Reagents
 Specimens
 Calibrators and Controls
 Performance Characteristics regarding
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
 - Method Comparison
 Reference Ranges
 Literature

- The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by
 DiaSys Diagnostic Systems GmbH & Co.KG
 Alte Strasse 9, 65558 Holzheim, Germany

#) Data entry by the user

*) Enter calibration or standard value and position

***) The factor must be checked by a calibration serum.