



MEGA CK-MB

PHOTOMETRIC TEST

Test ID Test name CK-MB

User code

Test code

				Report		
Normal range L/H	0/	25		Slope(%) Intercept A	100/	0
Panic range L/H	0/	2200		Slope(%) Intercept B	100/	0
Decimal point loc.	0			Unit	U/1	
Absorbance window	0/	0				

Calibrators

Calibration check	0 - 0	0.0000	Correction Clb.	N0.	0
Factor		8414,00			

Calibrator name	Conc.	Abs/Act.	(Factor)	S.V	Dil-Sample	Limits	L/H
C0 *****	0	0.000	()	10,0	0.0	0/	0
C1 _____	0	_____	()	0,0	0.0	0/	0
C2 _____	0	_____	()	0,0	0.0	0/	0
C3 _____	0	_____	()	0,0	0.0	0/	0
C4 _____	0	_____	()	0,0	0.0	0/	0
C5 _____	0	_____	()	0,0	0.0	0/	0
C6 _____	0	_____	()	0,0	0.0	0/	0
C7 _____	0	_____	()	0,0	0.0	0/	0
C8 _____	0	_____	()	0,0	0.0	0/	0

Test ID Test name CK-MB

User code

Test code

Assay

Measuring mode	2 points
Up/ Down	Up
Sample blank Test	CK-MB ()
Main/Sub W. L. (nm)	340/380
Calibration refer	_____ ()
Test read interval	76 -134
Blank read interval	0,0 - 0,0
Absorbance window	0,0/ 1,3
Max. Rate limit (%)	0

Sample

Determin. Per sample	1
Sample Vol./Dil. Vol.	10,0 0.0
Rerun/Dil. Vol	2,0 0.0

Regent

Regent	ID	Vol	H2O
Regent 1	_____	200	0
Regent 2	_____	50	0
Diluent	_____	0	0
Rerun Diluent	_____	0	0

Calibration

blank/Cal. - determin.	2/0
Type	Fator
Interval	0 Days 0 Hours

Reaction check

Endpoint check (Abs.)	0.000
Multiple determ. Range	0.000
Mono/Bi-chromatic	MAIN
Ratio/Differ check	Differ
End/Rate	End
Test read interval 1	0 - 0
Test read interval 2	0 - 0
Check limits L/H	0 - 0