

# MICROLAB 200

## TEST PARAMETERS

**CK-MB FS**

1	NAME	CK-MB
2	MODE	KINETIC
3	FILTER	340
4	TEMP	37°C
5	UNITS	U/L
6	R-BLANK	NO
	-L-Abs	0.000
	-H-Abs	2.300
1	CAL-NAME	
	.NR	
	.CONC	
	.REP	
	.L-Abs	0.000
	.H-Abs	2.300
2	FACTOR	8254**
3	S-BLANK	NO
1	DELAY	20
2	MEAS TIME	60
3	VOLUME	500
1	TEST REP	1
	.MAX DEV	0
2	REF-LOW	0
3	REF-HIGH	24
4	L-Abs	0.000
5	H-Abs	2.300
1	CONTROL-1	#
	.BATCH-NR	#
	.TARGET	#
	.LOW-LIM	#
	.HIGH-LIM	#
2	CONTROL-2	#
	.BATCH-NR	#
	.TARGET	#
	.LOW-LIM	#
	.HIGH-LIM	#

#) DATA ENTRY BY THE USER

\*) ENTER CALIBRATION OR STANDARD VALUE AND POSITION

\*\*) THE FACTOR MUST BE CHECKED BY A CALIBRATION SERUM.

### Order information

Cat. No.	Kit size
10 165 021	R1 5 x 20 ml + R2 1 x 25 ml
10 165 022	R1 5 x 80 ml + R2 1 x 100 ml
10 165 023	R1 1 x 800 ml + R2 1 x 200 ml

### Notes

1. Please refer to the package insert for CK-MB FS for the detailed information about the test on the following:

Clinical Relevance  
 Method and Principle  
 Composition and Stability of the Reagents  
 Specimens  
 Calibrators and Controls  
 Performance Characteristics concerning;  
   Measuring Range  
   Specificity/Interferences  
   Sensitivity/Limit of Detection  
   Precision (Reproducibility,  
   Repeatability)  
 Method Comparison  
 Reference Ranges  
 Literature

2. The stability of the reagent on board of the analyser is at least one month provided that contamination and evaporation are avoided.

3. Manufactured by  
 DiaSys Diagnostic Systems GmbH & Co.KG  
 Alte Strasse 9, 65558 Holzheim, Germany.