

BAYER OPERA

Resident Chemistry

CK-MB FS

NAME	CK
IMMUNOASSAY	No
TYPE	Zero Order
INVERSE CHEMISTRY	No
SAMPLE VOLUME	12.0
ALT. SAMPLE VOLUME	2.0
WAVELENGTH	340
DEPLETION TEST	Yes
BICHROMATIC CHEMISTRY	No
BICHROMATIC WAVELENGTH	-
BICHROMATIC TYPE	-
BICHROMATIC FACTOR 1	-
DEPLETION LIMIT	0.25
BICHROMATIC FACTOR 2	-
K1	-
K2	-
BICHROMATIC LIMIT 1	-
BICHROMATIC LIMIT 2	-
DELAY TIME	2:00
INCUBATION	-
BLANK TYPE	No
REAGENT VOLUME	320
SECOND REAGENT	80
2 ND REAGENT VOLUME	-
2 RGT DELAY	3:00
A1 DELAY	-
A2 DELAY	0:15
UNITS	U/l
UNIT FACTOR	1.000
DECIMAL POINT	0
RBL LOW	0.000
RBL HIGH	2:000
RANGE LOW	0
RANGE HIGH	1000
VALIDATION RANGE HIGH	1000
CALIBRATION FACTOR	**
REAGENT RATE	-
REAGENT BLANK	0.0
STANDARD VALUE	-
NORMAL LOW	0
NORMAL HIGH	24
SLOPE	1.000
INTERCEPT	0.000
ENDPOINT LIMIT	-
C1*10E-6	0.000
C2*10E-6	99999.01
D1*10E-6	90.0
DELTA	0.015
LINEARITY FACTOR	-
FIRST LIMIT	-
DAU METHOD	No
AUTO LINEARIZATION	No
CORRECTION LIMIT	-
AUTO LIN SLOPE	-
AUTO LIN INTERCEPT	-

TABLE IA	
No. OF STANDARDS	-
No. OF ASPIRATIONS	-
STANDARD 1	-
STANDARD 2	-
STANDARD 3	-
STANDARD 4	-
STANDARD 5	-
STANDARD 6	*

*) Enter calibrator or Standard value

#) Data entry by the user

**) Factor must be checked by using a calibrator

Order information

Cat. No.	Kit size
10 165 021	R1 5 x 20 ml + R2 1 x 25 ml
10 165 022	R1 5 x 80 ml + R2 1 x 100 ml
10 165 023	R1 1 x 800 ml + R2 1 x 200 ml

Notes

1. Please refer to the package insert for CK-MB FS for detailed information about the test on the following:

- Clinical Relevance
- Method and Principle
- Composition and Stability of the Reagents
- Specimens
- Calibrators and Controls
- Performance Characteristics regarding
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
 - Method Comparison
- Reference Ranges
- Literature

2. The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.

3. Manufactured by
DiaSys Diagnostic Systems GmbH
Alte Strasse 9, 65558 Holzheim, Germany