

BM / HITACHI 912

04-01 CHEMISTRY PARAMETERS

Enter the next code

CHLORIDE FS

TEST [CL]	[#]	TEST NAME [CL]	UNIT [mmol/l]
DATA MODE	[ON BOARD]	REPORT NAME [Chloride]	
CONTROL INT	[0]	INSTR. FACT. (Y=ax + b) a [1.0]	b [0.0]
EXPECTED VALUE	CLASS 1	EXPECTED VALUE	CLASS 2
AGE	M	F	
[] [] []	- [] []	- [] []	[] []
[] [] []	- [] []	- [] []	[] []
	[95] - [105]	[95] - [105]	[] []
TECHNICAL LIMIT	CLASS I	CLASS 2	
	[1] - [130]	[] - []	
STD	CON	POS	S.VOL PRE. DIL. VOL CODE LOT QUALITATIVE [NO]
(1)	[0.0][#]	[3]	[0] [0] [#] (1) [] []
(2)	[*] [#]	[3]	[0] [0] [#] (2) [] []
(3)	[] [] []	[] [] []	[] [] [] (3) [] []
(4)	[] [] []	[] [] []	[] [] [] (4) [] []
(5)	[] [] []	[] [] []	[] [] [] (5) [] []
(6)	[] [] []	[] [] []	[] [] [] (6) [] []

Order information

Cat. No.	Kit size
10 120 021	R 5 x 25 ml + 1 x 3 ml Std
1 1200 99 10 026	R 6 x 100 ml
10 120 023	R 1 x 1000 ml
10 120 0030	6 x 3 ml Standard
10 120 704	R 8 x 50 ml

Notes

- Please refer to the package insert for Chloride FS for the detailed information about the test on the following:

Clinical Relevance
 Method and Principle
 Composition and Stability of the Reagents
 Specimens
 Calibrators and Controls
 Performance Characteristics concerning:
 Measuring Range
 Specificity/Interferences
 Sensitivity/Limit of Detection
 Precision (Reproducibility, Repeatability)
 Method Comparison
 Reference Ranges
 Literature

- The stability of the reagent on board of the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by
 DiaSys Diagnostic Systems GmbH & Co.KG
 Alte Strasse 9, 65558 Holzheim, Germany.

TEST [CL]	
ASSAY CODE [1 point] [10] []	WAVELGTH (SUB / MAIN) [700] / [480]
ASSAY POINTS [31] - [0] - [0] - [0]	DILUTION [301][99] < CLASS 1 > < CLASS 2 >
S.VOL (NORMAL)	[3] [0] [0] [] [] []
S. VOL (DECREASE)	[2] [0] [0] [] [] []
S.VOL (INCREASE)	[6] [0] [0] [] [] []
ABS. LIMIT [0]	[] [2:INCREASE]
PROZONE LIMIT [32000]	[] [1:HIGHER]
REAGENT R1 [250] [0] [#] [#]	
R2 [0] [0] [#] [#]	
R3 [0] [0] [#] [#]	
R4 [0] [0] [#] [#]	
CALIB. TYPE [1:LINEAR]	[2] [2] [0] []
AUTO CALIB.	
TIME OUT BLANK [0]	SD LIMIT [0.1]
SPAN [0]	DUPLICATE LIMIT [100]
2 POINT [0]	SENSITIVITY LIMIT [0]
FULL [0]	SI ABS. LIMIT [-32000][32000]
CHANGE LOT [NO]	COMPENSATED LIMIT []
BOTTLE [NO]	

Data entry by the user
 * Enter calibration or standard value
 ** The given factor must be checked by a calibration serum.