

# CHLORIDE FS\*

## Order information

Cat. No.	Kit size
10 120 021	R 5 x 25 ml + 1 x 3 ml standard
1 1200 99 10 026	R 6 x 100 ml
10 120 023	R 1 x 1000 ml
10 120 030	6 x 3 ml standard

## Method

Colorimetric test, "Thiocyanate"

Chloride is releasing equivalent quantities of Thiocyanate from Mercury (II) Thiocyanate. Together with Iron ions Thiocyanate forms a red colored complex. The intensity of the color is proportional to the Chloride concentration.

## Reagent preparation and stability

The reagent is ready-to-use and stable up to the end of the indicated month of expiry, if contamination is avoided, stored at 15 – 25 °C and protected from light.

Stability of standard: until expiry date, when stored at 2 – 25 °C

## Specimen

Serum, plasma. Avoid hemolysis!

## Components and concentration in the test

Mercury (II) Thiocyanate	2 mmol/l
Mercury (II) Chloride	0.8 mmol/l
Ferric (III) Nitrate	20 mmol/l
Nitric acid	28 mmol/l

**Standard:** 100 mEq/l (mmol/l)

## Notes

The reagent contains Mercury (II) Thiocyanate and Mercury (II) Chloride. Do not swallow! Avoid contact with skin and mucous membranes.

## Normal range (see reference 1,2)

Adults:	98 - 110 mEq/l (98 - 110 mmol/l)
Children:	95 - 112 mEq/l (95 - 112 mmol/l)

## References

- Schoenfeld, R. G., Lewellen, C. J., Clin. Chem., 10, (1964), 533
- Witt, I., Trendelenburg, Chr., J. Clin. Chem. Clin. Biochem., 20, (1982), 235 - 242

\* fluid stable

# HITACHI 917

## Chemistry parameters

Test / Type	CL	Ser/PI
Assay / Time / Point	1 Point End A 10 A 20 0 0 0	
Wave (2 <sup>nd</sup> /Primary)	600 A 450 A	
S.Vol (Normal)	2 0 0	
S.Vol (Decrease)	1 0 0	
S.Vol (Increase)	4 0 0	
Diluent	00951 99	
Reagent (R1) T1	180 0 # 60	Timing R1
Reagent (R2) T2	0 0 # 0	R2
Reagent (R3) T3	0 0 # 0	R3
Reagent (R4) T4	0 0 # 0	
Abs. Limit	12000 Increase A	
Prozone Limit	32000 0 Lower A	
Cell Detergent	Detergent 1 A	

Calibration	
Calibration type	Linear A A
Point	2 Span Point 2
Weight	0
Autocalibration	
Blank	Time Out
Span	Change Over
2Point	Blank A
Full	Blank A
SD Limit	0.1
Duplicate limit	10 % 200 Abs
Sensitivity limit	-99999 999999
S1 Abs limit	-32000 32000

Range	
Application Code	# Unit mmol/l A
Report Name	Chloride
Data Mode	On Board A
Control Interval	1000
Instrument Factor (Y=aX+b)	a= 1.0 b= 0.0
Technical Limit	0 130
Repeat Limit	0 130
Expected Value	
(Male)	Qualitative
1 Y A	Cancel A
100 Y A	(1) 0
	(2) 0
	(3) 0
	(4) 0
	(5) 0
	(6) 0
(Female)	
1 Y A	
100 Y A	
(Default)	
Male A Range3 A	

Others	
<Standard>	(1) (2) (3) (4) (5) (6)
Calib. Code	501 #
Concentration	0 *
Position	
Sample Volume	2 2
Diluent S.Vol	0 0
Diluent Volume	0 0

#) Data entry by the user

\*) Enter calibration or standard value