

StarDust MC15

Chemistry parameters

CHLORIDE FS

Chloride

EP

| | |
|-------------|--------|
| WL1 | 500 |
| WL2 | 630 |
| Units | mmol/l |
| Temperature | 37 |
| Std | 100 |
| Factor | |
| Lin. Limit | 130 |
| Norm High | 108 |
| Norm Low | 96 |
| Decimals | 0 |
| Sample Vol | 5 |
| Reagent Vol | 500 |
| Reac. Time | 180 |

MC15

Procedure:

Reagent ready to use

Pipette into the cuvette strips in the appropriate way.

| | Reagent (Reading zone) | Sample (Sample zone) |
|----------------|---------------------------|-------------------------|
| Blank | 500 µl | |
| Std (Optional) | 500 µl | 5 µl |
| Samples | 500 µl | 5 µl |

Place the cuvette strip in the mixing position and press MIX

Place the cuvette strip in the incubating position for 3 minutes.

Go into programme no. 11 . Input the number of samples and the initial position for them and whether a recalibration is needed.

Place the cuvette strip in the reading position and press READ

Order information

| | | |
|------------------|----------|----------------------|
| Cat. No. | Kit size | |
| 10 120 021 | R 5 x | 25 ml + 1 x 3 ml Std |
| 1 1200 99 10 026 | R 6 x | 100 ml |
| 1 1200 99 10 730 | R 6 x | 20 ml |
| 10 120 030 | 6 x | 3 ml Standard |

Notes

1. Please refer to the package insert for Chloride FS for the detailed information about the test on the following:

- Clinical Relevance
- Method and Principle
- Composition and Stability of the Reagents
- Specimens
- Calibrators and Controls
- Performance Characteristics concerning:
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
- Method Comparison
- Reference Ranges
- Literature

2. Manufactured by
DiaSys Diagnostic Systems,
Alte Strasse 9, 65558 Holzheim, Germany.