

TEST PARAMETERS

Name	CHLORIDE		
Abbr. Name	CL		
Mode	ENDPOINT		
Wavelengths	506, 620		
Units	mmol/L		
Decimals	0		
Low Conc.	0 mmol/L		
High Conc.	130 mmol/L		
Calibrator Name	CAL		
Repeat	3		
Number	1		
Concentration	100 mmol/L		
Interval	7 days		
Cut off	NO		
Prozone Check	NO		
Ref. Male Low	98 mmol/L		
Ref Male High	107 mmol/L		
Ref Female Low	98		
Ref Female High	107		
Ref Ped. Low	98		
Ref Ped. High	107		
Control 1			
Control 2			
Control 3			
Correlat. Factor	1.000		
Correlat. Offset	0.000		
DUAL MODE	MONO MODE		
R1 bottle	25 ml	Sample blank	
Normal Volume	220µl	R1 bottle	25ml
Rerun Volume	220µl	Normal volume	220µl
Probe		Rerun volume	220µl
Sample Normal Vol.	3µl	Sample	
Sample Vol. rerun	1µl	Normal Volume	3µl
R2 bottle		rerun Volume	1µl
Normal Volume		Point one, two	
Rerun Volume		Incubation times	11.5 min
Predilution		Delay, min. time	
Incubation times	4.5 min	Linearity Limit	
Slope blank		Low Absorbance	0.00
Point one, two		High Absorbance	3.00
Delay, min. time		R.Abs L.Limit	0.00
Linearity Limit		R.Abs H. Limit	3.00
Low Absorbance	0.00	R.Abs Deviation	
High Absorbance	3.00	Reagent blank	Yes
R.Abs.L.Limit	0.00	Cal Low Limit	0
R.Abs. H.Limit	3.00	Cal. High Limit	0
R.Abs Diviation		Factor	1.00
Reagent blank	Yes		
Cal. Low Limit	0		
Cal.High Limit	0		
Factor	1.00		

Order information

Cat. No.	Kit size
10 120 021	R 5 x 25 ml + 1 x 3 ml Std

Notes

- Please refer to the package insert for Chloride FS for detailed information about the test on the following:

Clinical Relevance
 Method and Principle
 Composition and Stability of the Reagents
 Specimens
 Calibrators and Controls
 Performance Characteristics regarding

- Measuring Range
- Specificity/Interferences
- Sensitivity/Limit of Detection
- Precision (Reproducibility, Repeatability)
- Method Comparison

Reference Ranges
 Literature

- The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by
 DiaSys Diagnostic Systems GmbH & Co.KG
 Alte Strasse 9, 65558 Holzheim, Germany