

BM / HITACHI 912

04-01 CHEMISTRY PARAMETERS

CHOLESTEROL FS

TEST [CHOL] [#]	TEST NAME [CHOL]	UNIT [mg/dl]					
DATA MODE [ON BOARD]	REPORT NAME [CHOLESTEROL]						
CONTROL INT [0]	INSTR. FACT. (Y=ax + b) a [1.0]	b [0.0]					
EXPECTED VALUE	CLASS 1	EXPECTED VALUE CLASS 2					
AGE	M	F					
[] [] []	- [] []	- [] []					
[] [] []	- [] []	- [] []					
[0]	- [200]	[0] - [200]					
TECHNICAL LIMIT	CLASS I	CLASS 2					
[3]	- [750]	[] - []					
STD	CON	POS	S.VOL	PRE. DIL.	VOL	CODE LOT	QUALITATIVE
							[NO]
(1)	[0.0]	[#]	[3]	[0]	[0]	[#]	(1) [] []
(2)	[*]	[#]	[3]	[0]	[0]	[#]	(2) [] []
(3)	[]	[]	[]	[]	[]	[]	(3) [] []
(4)	[]	[]	[]	[]	[]	[]	(4) [] []
(5)	[]	[]	[]	[]	[]	[]	(5) [] []
(6)	[]	[]	[]	[]	[]	[]	(6) [] []

Order information

Cat. No.	Kit size
10 130 021	R 5 x 25 ml + 1 x 3 ml Std
1 1300 99 10 026	R 6 x 100 ml
10 130 023	R 1 x 1000 ml
10 130 030	6 x 3 ml Standard
10 130 704	R 8 x 50 ml

Notes

- Please refer to the package insert for Cholesterol FS for the detailed information about the test on the following:

Clinical Relevance
 Method and Principle
 Composition and Stability of the Reagents
 Specimens
 Calibrators and Controls
 Performance Characteristics concerning:
 Measuring Range
 Specificity/Interferences
 Sensitivity/Limit of Detection
 Precision (Reproducibility, Repeatability)
 Method Comparison
 Reference Ranges
 Literature

- The stability of the reagent on board of the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by
 DiaSys Diagnostic Systems GmbH & Co.KG
 Alte Strasse 9, 65558 Holzheim, Germany.

TEST [CHOL]
ASSAY CODE [1 point] [10] [] WAVELENGTH (SUB / MAIN) [700] / [505]
ASSAY POINTS [31] - [0] - [0] - [0] DILUTION [301] [99] < CLASS 1 > < CLASS 2 >
S.VOL (NORMAL) [3] [0] [0] [] [] []
S. VOL (DECREASE) [1] [0] [0] [] [] []
S.VOL (INCREASE) [6] [0] [0] [] [] []
ABS. LIMIT [0] [] [2:INCREASE]
PROZONE LIMIT [0] [] [1:HIGHER]
REAGENT R1 [250] [0] [#] [#]
R2 [0] [0] [#] [#]
R3 [0] [0] [#] [#]
R4 [0] [0] [#] [#]
CALIB. TYPE [1:LINEAR] [2] [2] [0] []
AUTO CALIB.
TIME OUT BLANK [0] SD LIMIT [0.1]
SPAN [0] DUPLICATE LIMIT [100]
2 POINT [0] SENSITIVITY LIMIT [0]
FULL [0] SI ABS. LIMIT [-32000][32000]
CHANGE LOT [NO] COMPENSATED LIMIT []
BOTTLE [NO]

Data entry by the user
 * Enter calibration or standard value
 ** The given factor must be checked by a calibration serum.
 ## Enter the next code