

StarDust MC15

Chemistry parameters

CHOLESTEROL FS

Chole

EP

WL1	500
WL2	630
Units	mg/dl
Temperature	37
Std	200
Factor	
Lin. Limit	380
Norm High	200
Norm Low	0
Decimals	0
Sample Vol	10
Reagent Vol	500
Reac. Time	600

MC15

Procedure:

Reagent ready to use

Pipette into the cuvette strips in the appropriate way.

	Reagent (Reading zone)	Sample (Sample zone)
Blank	500 µl	
Std (Optional)	500 µl	10 µl
Samples	500 µl	10 µl

Place the cuvette strip in the mixing position and press MIX

Place the cuvette strip in the incubating position for 10 minutes.

Go into programme no. 12 . Input the number of samples and the initial position for them and whether a recalibration is needed.

Place the cuvette strip in the reading position and press READ

Order information

Cat. No.	Kit size	
10 130 021	R 5 x	25 ml + 1 x 3 ml Std
1 1300 99 10 026	R 6 x	100 ml
10 130 023	R 1 x	1000 ml
10 130 030	6 x	3 ml Standard

Notes

1. Please refer to the package insert for Cholesterol FS for the detailed information about the test on the following

Clinical Relevance
Method and Principle
Composition and Stability of the Reagents
Specimens
Calibrators and Controls
Performance Characteristics concerning:
Measuring Range
Specificity/Interferences
Sensitivity/Limit of Detection
Precision (Reproducibility, Repeatability)
Method Comparison
Reference Ranges
Literature

2. Manufactured by
DiaSys Diagnostic Systems,
Alte Strasse 9, 65558 Holzheim, Germany