

BM / HITACHI 911

04-01 CHEMISTRY PARAMETERS

TEST [CHOL]	[#]	TEST NAME	[CHOL]	UNIT	[mg/dl]		
DATA MODE	[ON BOARD]	REPORT NAME	[CHOLESTEROL]				
CONTROL INT	[0]	INSTR. FACT.	(Y=aX + b)	a	[1.0]		
				b	[0.0]		
EXPECTED VALUE	CLASS 1	EXPECTED VALUE	CLASS 2				
AGE	M	F					
[]	[]	[]	- []	[]	[]		
[]	[]	[]	- []	[]	[]		
	[0]	- [200]	[0]	- [200]	[]		
TECHNICAL LIMIT	CLASS 1	CLASS 2					
	[3]	- [750]	[]	- []			
STD	CO	PO	S.VO	PRE. DIL.	VOL	CODE	QUALITATIVE
	N	S	L.			LOT	[NO]
(1)	[0.0]	[#]	[3]	[0]	[0]	[#]	(1) [] []
(2)	[*]	[#]	[3]	[0]	[0]	[#]	(2) [] []
(3)	[]	[]	[]	[]	[]	[]	(3) [] []
(4)	[]	[]	[]	[]	[]	[]	(4) [] []
(5)	[]	[]	[]	[]	[]	[]	(5) [] []
(6)	[]	[]	[]	[]	[]	[]	(6) [] []

TEST	[CHOL]				
ASSAY CODE	[1 point]	[10]	[]	WAVELGTH (SUB / MAIN)	[700] / [505]
ASSAY POINTS	[31] -	[0] -	[0] -	[0] DILUTION	[301] [99]
	< CLASS 1 >		< CLASS 2 >		
S.VOL (NORMAL)			[3]	[0]	[0] [] [] []
S. VOL (DECREASE)			[1]	[0]	[0] [] [] []
S.VOL (INCREASE)			[6]	[0]	[0] [] [] []
ABS. LIMIT	[0]		[]		[2:INCREASE]
PROZONE LIMIT	[0]		[]		[1:HIGHER]
REAGENT	R1	[250]	[0]	[#]	[#]
	R2	[0]	[0]	[#]	[#]
	R3	[0]	[0]	[#]	[#]
	R4	[0]	[0]	[#]	[#]
CALIB. TYPE		[1:LINEAR]		[2]	[2] [0] []
AUTOCALIB.					
TIME OUT BLANK	[0]		SD LIMIT		[0.1]
SPAN	[0]		DUPLICATE LIMIT		[100]
2 POINT	[0]		SENSITIVITY LIMIT		[0]
FULL	[0]		SI ABS. LIMIT		[-32000][32000]
CHANGE LOT	[NO]		COMPENSATED LIMIT		[]
BOTTLE	[NO]				

CHOLESTEROL FS

Order information

Cat. No.	Kit size	
10 130 021	R 5 x	25 ml + 1 x 3 ml Std
1 1300 99 10 026	R 6 x	100 ml
10 130 023	R 1 x	1000 ml
10 130 030	6 x	3 ml Standard
10 130 704	R 8 x	50 ml

Notes

- Please refer to the package insert for Cholesterol FS for detailed information about the test on the following:

Clinical Relevance
 Method and Principle
 Composition and Stability of the Reagents
 Specimens
 Calibrators and Controls
 Performance Characteristics regarding
 - Measuring Range
 - Specificity/Interferences
 - Sensitivity/Limit of Detection
 - Precision (Reproducibility, Repeatability)
 - Method Comparison
 Reference Ranges
 Literature

- The stability of the reagent on board the analyser is at least one month provided that contamination and evaporation are avoided.
- Manufactured by
 DiaSys Diagnostic Systems GmbH & Co.KG
 Alte Strasse 9, 65558 Holzheim, Germany

- # Data entry by the user
 * Enter calibration or standard value
 ** The given factor must be checked by a calibration serum.
 ## Enter the next code