

APPLICATION NOTE
SYNCHRON CX 4/5/7
COMPLEMENT C3 N-DIL (AUT-KIT)

1. Reagent Preparation

Sample: Ready for use

Buffer: Ready for use.

Antiserum: Ready for use.

Calibration: Dilute the Protein Standard High successively 1:2 in 9 g/L NaCl to set up a calibration curve. Use 9 g/L NaCl as zero point. Alternatively use the ready for use Protein Standard Set.

2. User Defined Chemistries

USER ID #	:				
Chemistry Name	:	C3			
Test Name	:	C3			
Reaction Type	:	Endpoint 2	Calculation Factor	:	0
Reaction Direction	:	Increase	Math. Model	:	1
Units	:	mg/dL	Cal. Time Limit	:	336 Hrs
Decimal Precision	:	X.X	No. of Calibrators	:	6
Primary Wavelength	:	340 nm	Secondary Wavelength	:	700 nm

Sample Volume	:	3 µL	Calibrators			
Primary inject Rgt	:		# 1	:	0.0	1 - 2 : 0.000
A	:	282 µL	# 2	:	*	2 - 3 : 0.000
			# 3	:	*	3 - 4 : 0.000
Secondary inject Rgt.	:		# 4	:	*	4 - 5 : 0.000
B	:	45 µL	# 5	:	*	5 - 6 : 0.000
Add Time	:	432 sec	# 6	:	*	6 - 1 : 0.000

REAGENT BLANK		REACTION	
Start Read	:	352 sec	Start Read : 680 sec
End Read	:	400 sec	End Read : 720 sec
Low ABS Limit	:	- 1.5	Low ABS Limit : - 1.5
High AS Limit	:	+1.5	High ABS Limit : +1.5

USABLE RANGE		SUBSTRATE DEPLETION	
Lower Limit	:	0.00	Initiate Rate : 99.9990
Upper Limit	:	300.00	Delta ABS : 1.5

RECOVERY / SENSITIVITY		
Std Dev (Conc.)	:	0.000
CV (%)	:	0.000
Std. Dev. (mA)	:	0.000
Threshold	:	0.000

* Standard Value

User defined

3. Order Information

C3C/AUT-000 1 x 10 mL Antiserum

5 x 25 mL Buffer

MPS/STH-001 Protein Standard High ,1 mL

MPS/STS-5X1 Protein Standard Set, 5 x 1 mL

MPC/CON-001 Protein Control, 1 mL

MPC/CON-005 Protein Control, 5 mL