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Amostra: Pronto para uso

Tampão R1: Pronto para uso

Anti-soro R2: Pronto para uso

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# ASLO (AUT – KIT)

### 1. Reagent preparation

Sample: Ready for use  
 Buffer: Ready for use  
 Latex: Ready for use  
 Standard: Dilute ASO Standard high successively 1:2 in saline 9 g/L to set up a calibration curve, use. Saline 9 g/L as zero point

### 2. Instrument settings

ASO			
General		Controls	
Mode	Differential – Bireagent	N° of controls	*
Unit	IU/mL	Replicates	*
Reaction type	Increasing	Rejection criteria	*
N° of replicates	1	Type	*
Decimals	1	Calculation model	*
<b>Procedure</b>		Name	*
Reading	Monochromatic	<b>Technical data</b>	
Principal	600	Lot	*
Secondary	-	Maximum value	*
Volumes		Minimum value	*
Sample	3	<b>Options</b>	
Reagent 1	280	Blanco with AD	No
Reagent 2	40	Blanc ABS limit	*
Wash	1.2	Linearity limit	*
Time		Detection limit	*
Measurement 1	Cycle N° 21 / 300 sec	Factor limit	*
Measurement 2	Cycle N° 41 / 600 sec	Reference range	*
Reagent 2	Cycle N° 22 / 315 sec		
Predilution factor	-		
Predilution type	-		
Post dilution factor	-		
Reduced	1		
Increased	1		
Automatic repetition	No		
<b>Calibration</b>			
Replicates blanco	1		
Replicates calibrator	1		
Calibration type	Specific		
Factor	-		
Concentration	**		
Calibration curve	Spline		
X	Logarithmic		
Y	Linear		

\* entered by operator

\*\* calibration concentrations

### 3. Ordering information

ASL/AUT-000            1 x 10 mL Antiserum  
                               5 x 25 mL Buffer  
 ASL/STH-001        ASLO standard high, 1 mL  
 ASL/CON-001        ASLO control, 1mL  
 139F003                Immunology Control Low, 1mL  
 139F002                Immunology Control High, 1 mL

# APPLICATION NOTE

## BIOSYSTEMS A25

# CRP (AUT – KIT)

### 1. Reagent preparation

Sample: Ready for use  
 Buffer: Ready for use  
 Antiserum: Ready for use  
 Standard: Dilute CRP Standard high successively 1:2 in saline 9 g/L to set up a calibration curve, use Saline 9 g/L as zero point.

### 2. Instrument settings

CRP			
General		Controls	
Mode	Differential – Bireagent	N° of controls	*
Unit	mg/dL	Replicates	*
Reaction type	Increasing	Rejection criteria	*
N° of replicates	1	Type	*
Decimals	1	Calculation model	*
<b>Procedure</b>		Name	*
Reading	Monochromatic	<b>Technical data</b>	
Principal	340	Lot	*
Secondary	-	Maximum value	*
Volumes		Minimum value	*
Sample	16	<b>Options</b>	
Reagent 1	250	Blanco with AD	No
Reagent 2	25	Blanc ABS limit	*
Wash	1.2	Linearity limit	*
Time		Detection limit	*
Measurement 1	Cycle N° 21 / 300 sec	Factor limit	*
Measurement 2	Cycle N° 41 / 600 sec	Reference range	*
Reagent 2	Cycle N° 22 / 315 sec		
Predilution factor	-		
Predilution type	-		
Post dilution factor	-		
Reduced	1		
Increased	1		
Automatic repetition	No		
<b>Calibration</b>			
Replicates blanco	1		
Replicates calibrator	1		
Calibration type	Specific		
Factor	-		
Concentration	**		
Calibration curve	Spline		
X	Logarithmic		
Y	Linear		

\* entered by operator

\*\* calibration concentrations

### 3. Ordering information

CRP/AUT-000      1 x 10 mL Antiserum  
                          5 x 25 mL Buffer  
 CRP/STH-001      CRP Standard High, 1 mL  
 CRP/COL-001      CRP Control Low, 1mL  
 CRP/COH-001      CRP Control High, 1 mL  
 139F003            Immunology Control Low, 1mL  
 139F002            Immunology Control High, 1 mL