

Quality Control

Certificate of Analysis

QC release date: 2018-04-05

Product: Bicarbonat FS

Product Code: 10950

Cat. No.: 1 0950 99 10 021

Delivery of kits: No of kits: 10
Kit Lot/Los: 60125417
Kit Expiry: 2019-05

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent Lot No.: 25155
Expiry: 2019-09

Examination:

Visual check: Reagent : clear, slightly yellow
Smell Reagent : no characteristic odour

pH: Reagent (25 °C ± 1 °C): 7.34

Absorbance: Reagent:
at 410 nm 1.224
at 415 nm 0.806

Control sera: TruLab Bicarbonat Lot.: 24830

Test results 19.9 mmol/L

Target value: 20.0 16.0 - 24.0

Additional remarks: Mf-Date: 21.03.2018

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka