

Quality Control

Certificate of Analysis

QC release date: 2017-12-08

Product: Bilirubin Auto Direct FS

Product Code: 10821

Cat. No.: 1 0821 99 10 021

Delivery of kits: No of kits: 9
Kit Lot/Los: 60124479
Kit Expiry: 2019-06

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent 1 Lot No.: 24518
Expiry: 2019-06

Reagent 2 Lot No.: 24518
Expiry: 2019-06

Examination:

Visual check Reagent 1: clear, colourless
Reagent 2: clear, colourless to yellow liquid

Smell Reagent 1: no characteristic odour
Reagent 2: no characteristic odour

pH: Reagent 1 (25 °C ± 1 °C): 1.15
Reagent 2 (25 °C ± 1 °C): 0.03

Absorbance: Reagent 1: 0.000
546 nm / 15-25 °C R1 + R2: 0.000

Control sera: TruLab N Trulab P
Lot: 22255 Lot: 22256

Test results 0.927 mg/dL 1.59 mg/dL

Target values 0.894 0.662 - 1.13 1.59 1.18 - 2.00

Additional remarks:

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka