

Quality Control

Certificate of Analysis

QC release date: 2019-01-11

Product: CK-NAC FS
Product Code: 11601

Cat. No.: 1 1601 99 10 021

Delivery of kits: No of kits: 6
Kit Lot/Los: 60130047
Kit Expiry: 2020-01

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent 1 Lot No.: 26664
Expiry: 2020-01

Reagent 2 Lot No.: 26664
Expiry: 2020-01

Examination:

Visual check Reagent 1: clear, colourless
Reagent 2: clear, colourless

Smell Reagent 1: slight smell (Thiol)
Reagent 2: no characteristic odour

pH: Reagent 1 (25 °C ± 1 °C): 6.05
Reagent 2 (25 °C ± 1 °C): 9.05
R1+R2 (37 °C ± 1 °C): 6.57

Absorbance: Reagent 1: 0.202
340 nm / 15-25 °C R1 + R2: 0.347

Control sera: TruLab N TruLab P
Lot: 24382 Lot: 24383

Fresh (at 37°C) 143 U/L 217 U/L

Target values 149 119 - 178 229 184 - 275

Additional remarks: Mf-Date: 14.01.2019

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed
by hand.

Sign.: Dr. G. Gorka