

Quality Control

## Certificate of Analysis

QC release date: 2018-04-27

**Product:** Calcium P FS  
**Product Code:** 11181

Cat. No.: 1 1181 99 10 920

---

**Delivery of kits:** No of kits: 45  
Kit Lot/Los: 60124896  
Kit Expiry: 2020-04

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent 1 Lot No.: 25156  
Expiry: 2020-04

Reagent 2 Lot No.: 25156  
Expiry: 2020-04

### Examination:

Visual check Reagent 1: clear, dark blue to violett

Reagent 2: clear, colourless

**pH:** Reagent 1 ( 25 °C ± 1 °C): 5.04  
Reagent 2 ( 25 °C ± 1 °C): 5.04

**Absorbance:** Reagent 1: 0.759  
660 nm / 37 °C Reagent 2: 0.000

**Control sera:** TruLab N TruLab P  
Lot: 22255 Lot: 22256

Test results 10.3 mg/dL 11.6 mg/dL

Target values 10.4 9.27 - 11.6 11.8 10.5 - 13.1

---

### Additional remarks:

This lot meets our specifications.

This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka