

Quality Control

Certificate of Analysis

QC release date: 2017-11-07

Product: Creatinine FS
Product Code: 11711

Cat. No.: 1 1711 99 10 920

Delivery of kits: No of kits: 15
Kit Lot/Los: 60121426
Kit Expiry: 2019-11

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent 1 Lot No.: 24520
Expiry: 2019-11

Reagent 2 Lot No.: 24520
Expiry: 2019-11

Examination:

Visual check Reagent 1: clear, colourless

Reagent 2: clear, yellow

Smell Reagent 1: no characteristic odour
Reagent 2: no characteristic odour

pH: Reagent 1 (25 °C ± 1 °C): 12.98
Reagent 2 (25 °C ± 1 °C): 1.74

Absorbance: Reagent 1: 0.000
492 nm / 25 °C R1 + R2: 0.176

Control sera: TruLab N TruLab P
Lot: 22255 Lot: 22256

Test results 1.59 mg/dL 4.30 mg/dL

Target values 1.60 1.25 - 1.95 4.32 3.37 - 5.26

Additional remarks:

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka