

Quality Control

## Certificate of Analysis

QC release date: 2019-04-16

**Product:** Creatinine PAP FS

**Product Code:** 11759

Cat. No.: 1 1759 99 10 026

---

**Delivery of kits:** No of kits: 2  
Kit Lot/Los: 60131812  
Kit Expiry: 2020-04

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent 1 Lot No.: 27077  
Expiry: 2020-04

Reagent 2 Lot No.: 27077  
Expiry: 2020-04

### Examination:

Visual check Reagent 1: clear, yellow  
Reagent 2: clear, pale brown

Smell Reagent 1: like detergent  
Reagent 2: no characteristic odour

**pH:** Reagent 1 ( 25 °C ± 1 °C): 8.09  
Reagent 2 ( 25 °C ± 1 °C): 8.03

**Absorbance:** Reagent 1: 0.004  
546 nm / 37 °C R1 + R2: 0.010

**Control sera:** TruLab N TruLab P  
Lot: 24382 Lot: 24383

Test results 1.50 mg/dL 4.90 mg/dL

Target values 1.46 1.14 - 1.78 4.88 3.81 - 5.95

---

**Additional remarks:** Mf-Date: 13.02.2019

This lot meets our specifications.  
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka