

Quality Control

Certificate of Analysis

QC release date: 2018-04-03

Product: Alkaline phosphatase FS (IFCC 37 °C)

Product Code: 10441

Cat. No.: 1 0441 99 10 021

Delivery of kits: No of kits: 74
Kit Lot/Los: 60126244
Kit Expiry: 2019-06

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent 1 Lot No.: 24932
Expiry: 2019-06

Reagent 2 Lot No.: 24932
Expiry: 2019-06

Examination:

Visual check Reagent 1: clear, colourless
Reagent 2: clear, pale yellow

Smell Reagent 1: slight smell (amines)
Reagent 2: no characteristic odour

pH: Reagent 1 (30 °C ± 1 °C): 10.41
Reagent 2 (25 °C ± 1 °C): 9.02
R1+R2 (30 °C ± 1 °C): 10.39

Absorbance: Reagent 1: 0.000
405 nm / 25 °C R1 + R2: 0.473

Control sera: TruLab N TruLab P
Lot: 22255 Lot: 22256

Fresh (at 37°C) 73.8 U/L 179 U/L

Target values 72.1 54.1 - 90.1 168 126 - 210

Additional remarks: Mf-Date: 21.03.2018

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed
by hand.

Sign.: Dr. G. Gorka