

Quality Control

## Certificate of Analysis

QC release date: 2017-11-24

**Product:** HDL-C Immuno FS  
**Product Code:** 13521

Cat. No.: 1 3521 99 10 026

---

**Delivery of kits:** No of kits: 7  
Kit Lot/Los: 60123951  
Kit Expiry: 2019-06

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent 1 Lot No.: 24888  
Expiry: 2019-06

Reagent 2 Lot No.: 24888  
Expiry: 2019-06

### Examination:

Visual check Reagent 1: clear, colourless to yellowish  
Reagent 2: clear, pale yellow

Smell Reagent 1: no characteristic odour  
Reagent 2: no characteristic odour

**pH:** Reagent 1 ( 25 °C ± 1 °C): 7.06  
Reagent 2 ( 25 °C ± 1 °C): 7.00

**Absorbance:** Reagent 1: 0.001  
600-700 nm / 37 °C R1 + R2: 0.003

**Control sera:** TruLab L 1 Trulab L 2  
Lot: 21238 Lot: 21239

Test results 49.8 mg/dL 43.4 mg/dL

Target values 50.4 40.3 - 60.5 43.5 34.8 - 52.2

---

### Additional remarks:

This lot meets our specifications.

This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka