

Quality Control

## Certificate of Analysis

QC release date: 2018-06-22

**Product:** LDH FS IFCC

**Product Code:** 14211

Cat. No.: 1 4211 99 10 920

---

**Delivery of kits:** No of kits: 2  
Kit Lot/Los: 60125806  
Kit Expiry: 2019-12

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent 1 Lot No.: 25431  
Expiry: 2019-12

Reagent 2 Lot No.: 25431  
Expiry: 2019-12

### Examination:

Visual check Reagent 1: clear, colourless  
Reagent 2: clear, yellowish

Smell Reagent 1: no characteristic odour  
Reagent 2: no characteristic odour

**pH:** Reagent 1 ( 37 °C ± 1 °C): 9.38  
Reagent 2 ( 25 °C ± 1 °C): 2.50  
R1+R2 ( 37 °C ± 1 °C): 9.32

**Absorbance:** Reagent 1: 0.153  
340 nm / 25 °C R1 + R2: 0.736

**Control sera:** TruLab N TruLab P  
Lot: 22255 Lot: 22256

Fresh (at 37°C) 157 U/L 252 U/L

Target values 157 129 - 185 262 215 - 309

---

**Additional remarks:** Mf-Date: 13.06.2018

This lot meets our specifications.  
This certificate of analysis is a computer printout and has therefore not been signed  
by hand.

Sign.: Dr. G. Gorka