

Quality Control

Certificate of Analysis

QC release date: 2018-09-10

Product: Total Protein UC FS
Product Code: 10210

Cat. No.: 1 0210 99 10 021

Delivery of kits: No of kits: 6
Kit Lot/Los: 60127562
Kit Expiry: 2020-04

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent Lot No.: 25699
Expiry: 2020-08

Examination:

Visual check: Reagent : clear, orange-red
Smell Reagent : no characteristic odour

pH: Reagent (25 °C ± 1 °C): 2.51

Absorbance: Reagent : 0.081
600 nm/ 15-25°C

Control sera:	TruLab Urine 1	TruLab Urine 2
	Lot: 23469	Lot: 23470
Fresh (at 37°C)	0.101 g/L	0.671 g/L
Target values	0.095 0.072 - 0.118	0.683 0.519 - 0.847

Additional remarks: Mf-Date: 15.08.2018

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka