

Quality Control

Certificate of Analysis

QC release date: 2020-02-03

Product: Total Protein UC FS
Product Code: 10210

Cat. No.: 1 0210 99 10 021

Delivery of kits: No of kits: 24
Kit Lot/Los: 60137622
Kit Expiry: 2021-08

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent Lot No.: 28286
Expiry: 2022-01-12

Examination:

Visual check: Reagent : clear, orange-red
Smell Reagent : no characteristic odour

pH: Reagent (25 °C ± 1 °C): 2.54

Absorbance: Reagent : 0.090
600 nm/ 15-25°C

Control sera:	TruLab Urine 1	TruLab Urine 2
	Lot: 26478	Lot: 26479
Fresh (at 37°C)	0.091 g/L	0.621 g/L
Target values	0.083 0.063 - 0.103	0.624 0.474 - 0.774

Additional remarks: Mf-Date: 15.01.2020

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka