

Quality Control

## Certificate of Analysis

QC release date: 2021-07-07

**Product:** Total Protein UC FS  
**Product Code:** 10210

Cat. No.: 1 0210 99 10 021

---

**Delivery of kits:** No of kits: 60  
Kit Lot/Los: 60146493  
Kit Expiry: 2023-02-26

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent Lot No.: 31062  
Expiry: 2023-06-11

### Examination:

Appearance: Reagent : clear, orange-red  
Odour: Reagent : no characteristic odour  
**pH:** Reagent ( 25 °C ± 1 °C): 2.53  
**Absorbance:** Reagent : 0.099  
600 nm/ 15-25°C

**Control sera:** TruLab Urine 1 TruLab Urine 2  
Lot: 30530 Lot: 30531  
Fresh (at 37°C) 0.079 g/L 0.668 g/L  
Target values 0.075 0.057 - 0.093 0.607 0.461 - 0.753

---

**Additional remarks:** Mf-Date: 11.06.2021

This lot meets our specifications.  
This certificate of analysis is a computer printout and has therefore not been signed  
by hand.

Sign.: Dr. G. Gorka