

Quality Control

## Certificate of Analysis

QC release date: 2020-09-25

**Product:** Triglycerides FS  
**Product Code:** 15710

Cat. No.: 1 5710 99 10 023

---

**Delivery of kits:** No of kits: 30  
Kit Lot/Los: 60141873  
Kit Expiry: 2022-09

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent Lot No.: 29206  
Expiry: 2022-09-04

### Examination:

Appearance: Reagent : clear, yellow

Odour: Reagent : like phenol

**pH:** Reagent ( 25 °C ± 1 °C): 7.17

**Absorbance:** Reagent : 0.020  
546 nm/ 15-25°C

Time to reach endpoint ≤ 10 minutes at 37 °C

<b>Control sera:</b>	TruLab N Lot: 25894	TruLab P Lot: 25895
Fresh (at 37°C)	94.9 mg/dL	179 mg/dL
Target values	92.9 76.2 - 110	171 141 - 202

---

**Additional remarks:** Mf-Date: 04.09.2020

This lot meets our specifications.  
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka