

Quality Control

Certificate of Analysis

QC release date: 2021-09-13

Product: Triglycerides FS
Product Code: 15710

Cat. No.: 1 5710 99 10 023

Delivery of kits: No of kits: 45
Kit Lot/Los: 60147823
Kit Expiry: 2023-09-03

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent Lot No.: 31069
Expiry: 2023-09-03

Examination:

Appearance: Reagent : clear, yellow

Odour: Reagent : like phenol

pH: Reagent (25 °C ± 1 °C): 7.19

Absorbance: Reagent : 0.016
546 nm/ 15-25°C

Time to reach endpoint ≤ 10 minutes at 37 °C

Control sera:	TruLab N Lot: 27547	TruLab P Lot: 27548
Fresh (at 37°C)	94.3 mg/dL	223 mg/dL
Target values	94.8 77.7 - 112	222 182 - 262

Additional remarks: Mf-Date: 03.09.2021

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka